	HTTRS, Inc		EPAIR FC D-PEN TON		
MODE	LS: Tono-Pen []	Tono-Pen 2 [] Tor	ıo-Pen XL [] Tono-P	en LTD [] Tono-Pen VET	[] Tono-Pen Avia []
c	OCU-LAB	BIO-RAD	MENTOR	MEDTRONIC	REICHERT
From:	Name: Address: Address: Phone:				
	Instrument(s) Serial N Problem(s): 275		HENDERSON NV	89074 • PHONE: (702) 8	396.9576

	TONOMETER REPAIR SERVICE		EPAIR FC D-PEN TON		and the second s				
MODELS: Tono-Pen [] Tono-Pen 2 [] Tono-Pen XL [] Tono-Pen LTD [] Tono-Pen VET [] Tono-Pen Avia []									
c	CU-LAB	BIO-RAD	MENTOR	MEDTRONIC	REICHERT				
From:	Name:								
	Address:								
	Phone:								
	Instrument(s) Serial N								
	Problem(s):								
	275	0 WOODBINE AVE	• HENDERSON NV	89074 • PHONE: (702) 89	96-9576				